

EPOS TAPE



IMPROVE YOUR PERFORMANCE



- About TS Co., Ltd Page 2
- What is Taping? Page 3
- EPOS ORIGINAL TAPE Page 4
- EPOS RAYON TAPE Page 5
- EPOS PRECUT TAPE Page 6, 7
- EPOS DESIGN TAPE Page 8
- EPOS TAPE SCISSORS Page 8
- EPOS CROSS TAPE Page 9
- Our Services Page 10
- Contact Page 11

ABOUT TS Co., Ltd

TS Co., Ltd

TS Co., Ltd. is a South Korean company based in Seoul and specialized in the production and distribution of kinesiology tapes. Our focus is to contribute to human's healthcare by serving the needs of our customers around the world with quality products.

Our activities:

Manufacturing and worldwide distribution of Kinesiology Taping Products

Our own brand is EPOS TAPE.

EPOS TAPE is manufactured in our own factory in Korea under ISO 9001 and ISO 13485 standards of Quality Management System, CE certification, Korea FDA registration and FDA registration in USA. Our taping products can be used in healthcare facilities, hospitals and pharmacies, as well as in the sport area.

Our R&D team constantly works on improving our offer in our own laboratory and our factory is equipped with the latest technology. It allows us to adapt and satisfy all the needs of our customers. Also, we have chosen only regional partners as suppliers of our raw materials in order to control all the production steps, from the supplies to the final products. This way we can ensure the best quality of our tapes.



WHAT IS TAPING?

By lifting the skin, EPOS KINESIOLOGY TAPE secures space for blood and lymph to flow smoothly again, allowing contusions and congestions to heal quickly.

Another important aspect of EPOS TAPE is the pain reduction. Our body can only accept one stimulus at a time and doesn't process further stimuli around the spinal cord once one has been taken into our body. Once applied, EPOS TAPE will constantly stimulate the sense of touch, which is processed faster than the sense of pain. Because of this, the sense of touch reaches the spinal cord first, blocking off the sense of pain. As a result, we feel less pain.


EPOS TAPE also gives your body more stability by supporting it with its elasticity equivalent to that of the human epidermis, providing an ideal grade of support. You can think of EPOS TAPE as a second, more resilient skin.

► Long Story Short

It lets you recover faster, can correct false and relieving postures and increases your performance by supporting joints and muscles. Also it can be used for cosmetic purposes such as curing bruises and wrinkles.

• PHYSIOLOGICAL EFFECTS OF EPOS TAPE •

EPOS TAPE features a special adhesive with a wave pattern design and an elasticity which is similar to the human skin. This offers a gentle and effective approach to the re-education of the neuromuscular system, the improvement of blood/lymph circulation and the relief of pain. Providing comfort and stability to your body.



1. EPOS ORIGINAL TAPE



- ▶ Main Functions
 - Normalizing muscle function
 - Improvement of lymphatic and blood flows
 - Pain reduction
 - Improvement of performance
 - Correction of the misalignment of joint

- ▶ Composition
 - Fabric: 97% Cotton, 3% Polyurethane
 - Adhesive: 100% Acrylic

- ▶ Features
 - Hypo-allergenic
 - Water Resistant
 - Breathable
 - Latex free
 - Elasticity similar as the one of the human skin and muscles
 - Last 3 days

- ▶ Origin - Made in South Korea
- ▶ Size - 5m x 5cm (2in x 16.4ft)
- ▶ Colors




2. EPOS RAYON TAPE



- ▶ Differences brought by the rayon fabric
 - a thin and shiny tape
 - a better elasticity
 - a smooth feeling when wearing
- ▶ Composition Fabric
 - 97% Rayon, 3% Polyurethane
 - Adhesive: 100% Acrylic
- ▶ Features
 - Hypo-allergenic
 - Water Resistant
 - Latex free
 - Elasticity similar as the one of the human skin and muscles
 - Last 3 days

- ▶ Origin - Made in South Korea
- ▶ Size - 5cm x 5m roll (2in x 16.4ft)
- ▶ Colors



3. EPOS PRECUT TAPE

• GENERAL INFORMATION •

• General information	<p>Designed and created to make the taping method accessible to individuals, EPOS PRECUT TAPE is convenient and easy to use.</p> <p>EPOS PRECUT TAPE contains:</p> <ul style="list-style-type: none">- High quality precut tape strips- Easy, fast and efficient applications- Detailed and clear illustrated instructions
• Composition	<p>Fabric: 97% Cotton, 3% Polyurethane Adhesive: 100% Acrylic</p>
• Features	<ul style="list-style-type: none">- Hypo-allergenic- Water Resistant- Latex free- Elasticity similar as the one of the human skin and muscles- Last 3 days
• Origin	Made in South Korea



EPOS TAPE

• EPOS PRECUT TAPE FOR THE UPPER BODY •



- Applications for the main painful zones of the upper body.
- Length of the precut strips adapted to the length of the muscles.

• 10 Precut strips



x2



x4



x4

X-shaped precut strips 20cm x 5cm-White I-shaped precut strips 20cm x 5cm-Red Y-shaped precut strips 25cm x 5cm-Black

• 12 Applications for the upper body

Neck, Shoulder, Golf Elbow, Tennis Elbow, Wrist, Biceps, Triceps, Upper Back 1, Upper Back 2, Lower Back, Pectoral, Abdomen

• EPOS PRECUT TAPE FOR THE LOWER BODY •



- Applications for the main painful zones of the lower body.
- Length of the precut strips adapted to the length of the muscles.

• 10 Precut strips



x3



x4



x3

Y-shaped precut strips 30cm x 5cm-White I-shaped precut strips 20cm x 5cm-Red Y-shaped precut strips 25cm x 5cm-Black

• 11 Applications for the lower body

Hip, Quadriceps, Inner Thigh, Hamstring, Knee Joint, Knee Joint Ligaments, Knee Joint Cartilage, Calf, Medial Ankle, Ankle, Achilles

4. EPOS DESIGN TAPE



✓ Blue Argyle

Combination of elegance with sport taping

✓ Lips

EPOS TAPE exclusive design

✓ EPOS Logo Tape

Help our distributors to promote EPOS TAPE


- ▶ Composition – Fabric: 97% Cotton, 3% Polyurethane
– Adhesive: 100% Acrylic
- ▶ Size – 5m x 5cm (2in x 16.4ft)

5. EPOS TAPE SCISSORS



- ▶ Features – Excellent sharpness
– Provides smooth cut
– No debris left on the blade after use
– Durable usage thanks to resilient materials
- ▶ Composition – Blade: Stainless steel with titanium coating
– Grip: Plastic
- ▶ Size – Total length: 21 cm – Blade length: 8 cm
- ▶ Package – Individual box with a protective sponge

6. EPOS CROSS TAPE




▶ Main Functions – EPOS CROSS TAPE can be used on acupuncture points and pain spots to increase the body self-healing powers

▶ Composition – Fabric: 100% Polyester
– Adhesive: 100% Acrylic

▶ Color – Beige


▶ Features – Hypo-allergenic and skin friendly
– Water resistant
– Latex free
– Resilient
– Non elastic
– Last 3 days

size A




- 21mm x 27mm
- 9 patches per sheet

size B



- 28mm x 36mm
- 6 patches per sheet

size C



- 44mm x 52mm
- 2 patches per sheet

OUR *SERVICES*

1. *OEM*

We can produce your own brand of Kinesiology Tape or Cooling Gel Sheet. Contact us for more information.

2. *Size of Kinesiology taping rolls*

We can adapt the size of our rolls.



◎ CONTACT

✓ Websites

www.epostape.net

✓ Address

#900 Industrial Complex, Hagun-Ri, Yangchon-Eup,
Gimpo-Si, Gyeonggi-Do 415-843, South Korea

✓ Phone number

+82 70 8768 7064


✓ E-mail

master@epostape.net

✓ EPOS TAPE Community

 www.facebook.com/epostape

 [@epostape](https://twitter.com/epostape)

 epos-tape.blogspot.kr/

 www.flickr.com/photos/epos-tape/



E EPOS TAPE